

Product datasheet for **TA331719**

Nucleotide binding protein like (NUBPL) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-NUBPL Antibody is: synthetic peptide directed towards the C-terminal region of Human NUBPL. Synthetic peptide located within the following region: AQTLGLEVLGDIPLHLNIREASDTGQPIVFSQPESDEAKAYLRIAVEVVR
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	32 kDa
Gene Name:	nucleotide binding protein like
Database Link:	NP_079428 Entrez Gene 80224 Human Q8TB37
Background:	This gene encodes a member of the Mrp/NBP35 ATP-binding proteins family. The encoded protein is required for the assembly of the respiratory chain NADH dehydrogenase (complex I), an oligomeric enzymatic complex located in the inner mitochondrial membrane. The respiratory complex I consists of 45 subunits and 8 iron-sulfur (Fe/S) clusters. This protein is an Fe/S protein that plays a critical role in the assembly of respiratory complex I, likely by transferring Fe/S into the Fe/S-containing complex I subunits. Mutations in this gene cause mitochondrial complex I deficiency. Alternatively spliced transcript variants encoding distinct isoforms have been identified.



[View online »](#)

Synonyms: C14orf127; hulnd1; IND1

Note: Immunogen sequence homology: Dog: 100%; Rat: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Pig: 93%; Guinea pig: 93%; Mouse: 92%; Rabbit: 92%

Product images:



Host: Rabbit; Target Name: NUBPL; Sample Tissue: A549 Whole Cell lysates; Antibody Dilution: 1.0ug/ml; NUBPL is supported by BioGPS gene expression data to be expressed in A549