

Product datasheet for **TA331753**

HROB Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-C17orf53 Antibody is: synthetic peptide directed towards the middle region of Human C17orf53. Synthetic peptide located within the following region: QPFQSPSSWLSGKAHLPRPRTPNSSCSTPSRTSSGLFPRIPLQPQAPVSS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	57 kDa
Gene Name:	chromosome 17 open reading frame 53
Database Link:	NP_076937 Entrez Gene 78995 Human Q8N3J3
Background:	The function of this protein remains unknown.
Synonyms:	FLJ11594; MGC3130
Note:	Immunogen sequence homology: Human: 100%; Horse: 92%; Dog: 86%; Rabbit: 86%; Bovine: 79%



[View online »](#)

Product images:



Host: Rabbit; Target Name: C17orf53; Sample
Tissue: OVCAR-3 Whole Cell lysates; Antibody
Dilution: 1.0ug/ml