

Product datasheet for **TA331847**

ZNF313 (RNF114) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ZNF313 Antibody: synthetic peptide directed towards the N terminal of human ZNF313. Synthetic peptide located within the following region: AVELERQIESTETSCHGCRKNFFLSKIRSHVATCSKYQNYIMEGVKATIK
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	26 kDa
Gene Name:	ring finger protein 114
Database Link:	NP_061153 Entrez Gene 55905 Human Q9Y508
Background:	ZNF313 may play a role in spermatogenesis
Synonyms:	PSORS12; ZNF313
Note:	Immunogen sequence homology: Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Rabbit: 100%; Dog: 93%; Mouse: 93%; Guinea pig: 92%; Yeast: 75%



[View online »](#)

Product images:



WB Suggested Anti-ZNF313 Antibody Titration:
1.25ug/ml; ELISA Titer: 1: 312500; Positive
Control: Transfected 293T