

Product datasheet for **TA331874**

ZNF312 (FEZF2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ZNF312 Antibody: synthetic peptide directed towards the N terminal of human ZNF312. Synthetic peptide located within the following region: MASSASLETMVPPACPRAGASPATSKTLAFSIERIMAKTSEPRAPFPRP
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	50 kDa
Gene Name:	FEZ family zinc finger 2
Database Link:	NP_060478 Entrez Gene 55079 Human Q8TBJ5
Background:	ZNF312 is a candidate transcription factor
Synonyms:	FEZ; FEZL; FKSG36; TOF; ZFP312; ZNF312
Note:	Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%



[View online »](#)

Product images:



WB Suggested Anti-ZNF312 Antibody Titration:
2.5ug/ml; Positive Control: Jurkat cell lysate