

Product datasheet for **TA332055**

ABCF2 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ABCF2 Antibody: synthetic peptide directed towards the N terminal of human ABCF2. Synthetic peptide located within the following region: DIYHLTREMPSPDKTPLHCVMEVDTERAMLEKEAERLAHEDAEECKLMEL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	72 kDa
Gene Name:	ATP binding cassette subfamily F member 2
Database Link:	NP_005683 Entrez Gene 114483834 Human Entrez Gene 10061 Human Q9UG63
Background:	This gene encodes a member of the ATP-binding cassette (ABC) transporter superfamily. ATP-binding cassette proteins transport various molecules across extra- and intracellular membranes. Alterations in this gene may be involved in cancer progression. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 3 and 7. [provided by RefSeq, Jul 2013]
Synonyms:	ABC28; EST133090; HUSSY-18; HUSSY18



[View online »](#)

Note: Immunogen sequence homology: Pig: 100%; Human: 100%; Guinea pig: 100%; Dog: 93%; Rat: 93%; Horse: 93%; Mouse: 93%; Bovine: 93%; Rabbit: 93%; Zebrafish: 93%

Product images:



WB Suggested Anti-ABCF2 Antibody Titration: 0.2-1 ug/ml; Positive Control: Human kidney