

Product datasheet for **TA332177**

BTBD18 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-BTBD18 Antibody is: synthetic peptide directed towards the C-terminal region of Human BTBD18. Synthetic peptide located within the following region: IEGEEWCLPDMELWPRELTELEKEPAGENRGPTSELLSPLVMPSEVSEVLS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	78 kDa
Gene Name:	BTB domain containing 18
Database Link:	NP_001138573 Entrez Gene 643376 Human B2RXH4
Background:	The function of this protein remains unknown.
Synonyms:	BTBD18
Note:	Immunogen sequence homology: Human: 100%; Pig: 79%



[View online »](#)

Product images:



Host: Rabbit; Target Name: BTBD18; Sample Tissue: 721_B Whole Cell lysates; Antibody Dilution: 1.0ug/ml; BTBD18 is supported by BioGPS gene expression data to be expressed in 721_B