

Product datasheet for **TA332342**

OR8K1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-OR8K1 Antibody is: synthetic peptide directed towards the C-terminal region of Human OR8K1. Synthetic peptide located within the following region: ILRMNSRKGRYKAFSTCSSHLTVWIMFYGTLLFIYLQPKSSHTLAIDKMA
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	35 kDa
Gene Name:	olfactory receptor family 8 subfamily K member 1
Database Link:	NP_001002907 Entrez Gene 390157 Human Q8NGG5
Background:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms.



[View online »](#)

Synonyms:	OR8N1P; OR11-182
Note:	Immunogen sequence homology: Human: 100%; Bovine: 93%; Pig: 91%; Rat: 91%; Mouse: 91%; Dog: 86%; Horse: 86%; Rabbit: 86%; Guinea pig: 85%
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction

Product images:

Host: Rabbit
Target Name: OR8K1
Sample Tissue: HepG2 Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit; Target Name: OR8K1; Sample Tissue: HepG2 Whole Cell lysates; Antibody Dilution: 1.0ug/ml