

Product datasheet for **TA332350**


IDNK Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-IDNK Antibody is: synthetic peptide directed towards the middle region of Human IDNK. Synthetic peptide located within the following region: RVVLACSALKKTYRDILTQGKDGVALKCEESGKEAKQAEMQLLVHLSGS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	25 kDa
Gene Name:	idnK, gluconokinase homolog (E. coli)
Database Link:	NP_001001551 Entrez Gene 414328 Human Q5T6J7
Background:	The function of this protein remains unknown.
Synonyms:	bA522I20.2; C9orf103; hGntK
Note:	Immunogen sequence homology: Human: 100%



[View online »](#)

Product images:

Host: Rabbit; Target Name: IDNK; Sample Tissue: MCF7 Whole Cell lysates; Antibody Dilution: 1.0ug/ml