

Product datasheet for **TA332355**

Monoamine Oxidase B (MAOB) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-MAOB Antibody: synthetic peptide directed towards the C terminal of human MAOB. Synthetic peptide located within the following region: GKIPEDIWQSEPESVDVPAQPITTTFLERHLPSVPGLLRLIGLTTIFSA
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	59 kDa
Gene Name:	monoamine oxidase B
Database Link:	NP_000889 Entrez Gene 4129 Human P27338



[View online »](#)

Background: MAOB belongs to the flavin monoamine oxidase family. It is a enzyme located in the mitochondrial outer membrane. It catalyzes the oxidative deamination of biogenic and xenobiotic amines and plays an important role in the metabolism of neuroactive and vasoactive amines in the central nervous ssystem and peripheral tissues. This protein preferentially degrades benzylamine and phenylethylamine. The protein encoded by this gene belongs to the flavin monoamine oxidase family. It is a enzyme located in the mitochondrial outer membrane. It catalyzes the oxidative deamination of biogenic and xenobiotic amines and plays an important role in the metabolism of neuroactive and vasoactive amines in the central nervous ssystem and peripheral tissues. This protein preferentially degrades benzylamine and phenylethylamine.


Synonyms: MGC26382

Note: Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Guinea pig: 100%; Rabbit: 93%; Bovine: 92%

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Arginine and proline metabolism, Drug metabolism - cytochrome P450, Glycine, serine and threonine metabolism, Histidine metabolism, Metabolic pathways, Phenylalanine metabolism, Tryptophan metabolism, Tyrosine metabolism

Product images:



WB Suggested Anti-MAOB Antibody Titration: 0.2-1 ug/ml; Positive Control: Human Liver