

Product datasheet for TA332358

TAP1 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: The immunogen for Anti-TAP1 Antibody: synthetic peptide directed towards the middle

region of human TAP1. Synthetic peptide located within the following region:

LVTFVLYQMQFTQAVEVLLSIYPRVQKAVGSSEKIFEYLDRTPRCPPSGL

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 87 kDa

Gene Name: transporter 1, ATP-binding cassette, sub-family B (MDR/TAP)

Database Link: NP 000584

Entrez Gene 6890 Human

Q03518



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Background:

The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). TAP1 is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance. TAP1 is involved in the pumping of degraded cytosolic peptides across the endoplasmic reticulum into the membrane-bound compartment where class I molecules assemble. The membraneassociated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance. The protein encoded by this gene is involved in the pumping of degraded cytosolic peptides across the endoplasmic reticulum into the membrane-bound compartment where class I molecules assemble. Mutations in this gene may be associated with ankylosing spondylitis, insulin-dependent diabetes mellitus, and celiac disease. Publication Note: This RefSeg record includes a subset of the publications that are available for this gene. Please see the Entrez Gene record to access additional publications.

Synonyms: ABC17; ABCB2; APT1; D6S114E; PSF-1; PSF1; RING4; TAP1 0102N; TAP1N

Note: Immunogen sequence homology: Rat: 100%; Human: 100%; Mouse: 100%; Horse: 93%;

Bovine: 92%; Rabbit: 92%; Pig: 86%; Guinea pig: 86%; Zebrafish: 83%

Protein Families: Druggable Genome, Transmembrane


Protein Pathways: ABC transporters, Antigen processing and presentation, Primary immunodeficiency

Product images:





WB Suggested Anti-TAP1 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive Control: MCF7 cell lysate






Host: Rabbit; Target Name: TAP1; Sample Tissue: Hela; Antibody Dilution: 1.0ug/ml; TAP1is supported by BioGPS gene expression data to be expressed in HeLa







Rabbit Anti-TAP1 Antibody; Formalin Fixed Paraffin Embedded Tissue: Human Lung Tissue; Observed Staining: Cytoplasmic in alveolar type I & II cells; Primary Antibody Concentration: 1: 100; Secondary Antibody: Donkey anti-Rabbit-Cy3; Secondary Antibody Con