

## Product datasheet for **TA332418**

### PDPK1 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB 1:500 - 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Protein of human PDPK1
Formulation:	Store at -20°C (regular) and -80°C (long term). Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	3-phosphoinositide dependent protein kinase 1
Database Link:	<a href="#">NP_002604</a> <a href="#">Entrez Gene 18607 Mouse</a> <a href="#">Entrez Gene 81745 Rat</a> <a href="#">Entrez Gene 5170 Human</a> <a href="#">O15530</a>
Synonyms:	PDK1; PDPK2; PDPK2P; PRO0461
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Endometrial cancer, Focal adhesion, Insulin signaling pathway, mTOR signaling pathway, Non-small cell lung cancer, PPAR signaling pathway, Prostate cancer



[View online »](#)

## Product images:



Western blot analysis of extracts of various cell lines, using PDPK1 antibody.