

Product datasheet for **TA332485S**

HSD17B2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB 1:500 - 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human HSD17B2
Formulation:	Store at -20°C (regular) and -80°C (long term). Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	42 kDa
Gene Name:	hydroxysteroid (17-beta) dehydrogenase 2
Database Link:	NP_002144 Entrez Gene 15486 Mouse Entrez Gene 3294 Human P37059
Synonyms:	EDH17B2; HSD17; SDR9C2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Androgen and estrogen metabolism, Metabolic pathways



[View online »](#)

Product images:



Western blot analysis of extracts of various cell lines, using HSD17B2 antibody.