

Product datasheet for **TA332689S**

MERTK Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ICC/IF, IHC, WB
Recommended Dilution:	WB 1:500 - 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human MERTK
Formulation:	Store at -20°C (regular) and -80°C (long term). Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	999
Gene Name:	MER proto-oncogene, tyrosine kinase
Database Link:	NP_006334 Entrez Gene 17289 MouseEntrez Gene 65037 RatEntrez Gene 10461 Human Q12866
Synonyms:	c-mer; MER; RP38
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane



[View online »](#)

Product images:



Western blot analysis of extracts of various cell lines, using MERTK antibody.