

Product datasheet for **TA332775**


Rabl2 Guinea Pig Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB 1:500 - 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Guinea Pig
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human Rabl2
Formulation:	Store at -20°C (regular) and -80°C (long term). Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	228
Gene Name:	RAB, member RAS oncogene family-like 2
Database Link:	NP_081093 Entrez Gene 362987 Rat Entrez Gene 68708 Mouse E9Q9D5
Synonyms:	1110031N17Rik; AW549303; Rabl2a



[View online »](#)

Product images:

Western blot analysis of recombinant protein using Rab12 antibody.