

## Product datasheet for **TA332778**

### Eif2ak2 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB 1:500 - 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human Eif2ak2
Formulation:	Store at -20°C (regular) and -80°C (long term). Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	515
Gene Name:	eukaryotic translation initiation factor 2-alpha kinase 2
Database Link:	<a href="#">NP_035293</a> <a href="#">Entrez Gene 5610 Human</a> <a href="#">Entrez Gene 54287 Rat</a> <a href="#">Entrez Gene 19106 Mouse</a> <a href="#">Q03963</a>
Synonyms:	EIF2AK1; MGC126524; OTTHUMP00000201320; PKR; PRKR



[View online »](#)

## Product images:



Western blot analysis of extracts of MCF7 cell line, using Eif2ak2 antibody.