

Product datasheet for TA332891

Product datasileet for TASS2691

53BP1 (TP53BP1) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB 1:500 - 1:2000

Reactivity: Human, Mouse, Rat

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: Recombinant protein of human RAD51

Formulation: Store at -20°C (regular) and -80°C (long term). Avoid freeze / thaw cycles. Buffer: PBS with

0.02% sodium azide, 50% glycerol, pH7.3.

Concentration: lot specific

Purification: Affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: tumor protein p53 binding protein 1

Database Link: NP 005648

Entrez Gene 296099 RatEntrez Gene 7158 Human

Q12888

Synonyms: 53BP1; p202

Protein Families: Druggable Genome, Transcription Factors




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Product images:



Western blot analysis of extracts of A431 cell line, using TP53BP1 antibody.