

Product datasheet for **TA332920**

Pea3 (ETV4) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ICC/IF, IHC, WB
Recommended Dilution:	WB 1:500 - 1:2000;IF 1:50- 1:200
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human ETV4
Formulation:	Store at -20°C (regular) and -80°C (long term). Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	484
Gene Name:	ETS variant 4
Database Link:	NP_001977 Entrez Gene 18612 MouseEntrez Gene 360635 RatEntrez Gene 2118 Human P43268
Synonyms:	E1A-F; E1AF; PEA3; PEAS3
Protein Families:	Druggable Genome, Transcription Factors




[View online »](#)

Product images:



Western blot analysis of extracts of various cell lines, using ETV4 antibody.



Immunofluorescence analysis of U2OS cell using ETV4 antibody. Blue: DAPI for nuclear staining.