

Product datasheet for **TA332980**

VASH1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB 1:500 - 1:2000
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human VASH1
Formulation:	Store at -20°C (regular) and -80°C (long term). Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	37 kDa
Gene Name:	vasohibin 1
Database Link:	NP_055724 Entrez Gene 238328 Mouse Entrez Gene 22846 Human Q7L8A9
Synonyms:	KIAA1036
Protein Families:	Secreted Protein



[View online »](#)

Product images:



Western blot analysis of extracts of various cell lines, using VASH1 antibody.