

Product datasheet for **TA333321**


EZHIP Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-CXorf67 Antibody: synthetic peptide directed towards the middle region of human LOC340602. Synthetic peptide located within the following region: RRSLSGSADENPSCGTGSRERLAFQSRSGSPDPEVPSRASPPVWHAVRMRA
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	52 kDa
Gene Name:	chromosome X open reading frame 67
Database Link:	NP_981952 Entrez Gene 340602 Human Q86X51
Background:	The specific function of the protein remains unknown.
Synonyms:	MGC47837
Note:	Immunogen sequence homology: Human: 100%; Dog: 77%



[View online »](#)

Product images:

WB Suggested Anti-CXorf67 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive
Control: Transfected 293T