

Product datasheet for **TA333372**

KANK3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ANKRD47 Antibody: synthetic peptide directed towards the N terminal of human ANKRD47. Synthetic peptide located within the following region: MAKFALNQNLPDLGGPRLCPVPAAGGARSPSPYSVETPYGFHLDLDFLK
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	86 kDa
Gene Name:	KN motif and ankyrin repeat domains 3
Database Link:	NP_940873 Entrez Gene 256949 Human Q6NY19
Background:	The specific function of this protein remains unknown.
Synonyms:	ANKRD47
Note:	Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%



[View online »](#)

Product images:



WB Suggested Anti-ANKRD47 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive
Control: Human heart