

Product datasheet for TA333387


WDR53 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-WDR53 Antibody: synthetic peptide directed towards the middle region of human WDR53. Synthetic peptide located within the following region: NLLASADDSGAIKILDLENKKVIRSLKRHSNICSSVAFRPQRPQSLVSCG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	39 kDa
Gene Name:	WD repeat domain 53
Database Link:	NP_872433 Entrez Gene 348793 Human Q7Z5U6
Background:	The specific function of this protein remains unknown.
Synonyms:	MGC12928; MGC64882
Note:	Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Guinea pig: 100%; Rabbit: 92%; Zebrafish: 85%


[View online »](#)

Product images:



WB Suggested Anti-WDR53 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:1562500; Positive
Control: Human brain