

Product datasheet for **TA333426**

C2orf53 (PRR30) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-C2orf53 Antibody: synthetic peptide directed towards the middle region of human C2orf53. Synthetic peptide located within the following region: PQGKATQACGHQLPASQPPAAQARADPVPGPSQTRSFERSAGLQSPNSPR
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	45 kDa
Gene Name:	proline rich 30
Database Link:	NP_848648 Entrez Gene 339779 Human Q53SZ7
Background:	The exact functions of C2orf53 remain unknown.
Synonyms:	C2orf53
Note:	Immunogen sequence homology: Human: 100%



[View online »](#)

Product images:



WB Suggested Anti-C2orf53 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:1562500; Positive
Control: 721_B cell lysate