

Product datasheet for **TA333486**


C1orf172 (KDF1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-C1orf172 Antibody is: synthetic peptide directed towards the C-terminal region of Human C1orf172. Synthetic peptide located within the following region: LSQDELTVQISQETTADAIARKLRPYGAPGYASHDSSFQGTDTDSSGAP
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	43 kDa
Gene Name:	keratinocyte differentiation factor 1
Database Link:	NP_689578 Entrez Gene 126695 Human Q8NAX2
Background:	The function of this protein remains unknown.
Synonyms:	C1orf172
Note:	Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Mouse: 93%



[View online »](#)

Product images:

Host: Rabbit; Target Name: C1orf172; Sample Tissue: Fetal Brain lysates; Antibody Dilution: 1.0 ug/ml.