

Product datasheet for **TA333782**

C1orf41 (HSPB11) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-HSPB11 Antibody is: synthetic peptide directed towards the middle region of Human HSPB11. Synthetic peptide located within the following region: IERLVIQSYFVQTLKIEKSTSKEPVDFEQWIEKDLVHTEGQLQNEEIVAH
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	13 kDa
Gene Name:	heat shock protein family B (small) member 11
Database Link:	NP_057210 Entrez Gene 51668 Human Q9Y547
Background:	The function of this protein remains unknown.
Synonyms:	C1orf41; HSPCO34; IFT25; PP25
Note:	Immunogen sequence homology: Dog: 100%; Human: 100%; Horse: 93%; Bovine: 93%; Rabbit: 93%; Pig: 86%; Rat: 85%; Mouse: 79%



[View online »](#)

Product images:

WB Suggested Anti-HSPB11 Antibody; Titration:
1.0 ug/ml; Positive Control: ACHN Whole Cell