

Product datasheet for **TA333857**

MAPKBP1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-MAPKBP1 Antibody: synthetic peptide directed towards the N terminal of human MAPKBP1. Synthetic peptide located within the following region: MAVEGSTITSRIKNLLRSPSIKLRRSKAGNRREDLSSKVTLEKVLGITVS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	163 kDa
Gene Name:	mitogen-activated protein kinase binding protein 1
Database Link:	NP_055809 Entrez Gene 23005 Human O60336
Background:	MAPKBP1 is involved in JNK signaling pathway.
Synonyms:	JNKBP-1
Note:	Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%
Protein Families:	Druggable Genome



[View online »](#)

Product images:



WB Suggested Anti-MAPKBP1 Antibody Titration:
0.2-1 ug/ml; Positive Control: Human Stomach