

Product datasheet for **TA333999**

ZNF202 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ZNF202 Antibody: synthetic peptide directed towards the middle region of human ZNF202. Synthetic peptide located within the following region: LDPTQKEFYGEYVLEEDCGIVVLSLSPFIPRPDEISQVREEEPWVPDIQEP
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	75 kDa
Gene Name:	zinc finger protein 202
Database Link:	NP_003446 Entrez Gene 7753 Human O95125



[View online »](#)


Background: ZNF202 is a transcriptional repressor that binds to elements found predominantly in genes that participate in lipid metabolism. Among its targets are structural components of lipoprotein particles (apolipoproteins AIV, CIII, and E), enzymes involved in lipid processing (lipoprotein lipase, lecithin cholesteryl ester transferase), transporters involved in lipid homeostasis (ABCA1, ABCG1), and several genes involved in processes related to energy metabolism and vascular disease. Western blots using two different antibodies against two unique regions of this protein target confirm the same apparent molecular weight in our tests.

Synonyms: ZKSCAN10; ZSCAN42

Note: Immunogen sequence homology: Human: 100%; Dog: 93%; Horse: 93%; Bovine: 93%; Pig: 92%; Rat: 86%; Guinea pig: 86%; Rabbit: 85%; Mouse: 79%

Protein Families: Transcription Factors

Product images:



WB Suggested Anti-ZNF202 Antibody Titration: 0.2-1 ug/ml; Positive Control: Jurkat cell lysate. ZNF202 is supported by BioGPS gene expression data to be expressed in Jurkat