

## Product datasheet for **TA334183**


### DMRT2 Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 2-4ug/mL
Reactivity:	Human (Expected from sequence similarity: Mouse, Rat, Cow)
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Peptide with sequence C-KRPSSAITRVSQ, from the C Terminus of the protein sequence according to NP_006548.1; NP_870987.2.
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	doublesex and mab-3 related transcription factor 2
Database Link:	<a href="#">NP_870987</a> <a href="#">Entrez Gene 226049 MouseEntrez Gene 309430 RatEntrez Gene 10655 Human Q9Y5R5</a>
Synonyms:	doublesex and mab-3 related transcription factor 2; OTTHUMP00000020964; OTTHUMP00000020965; OTTHUMP00000020966; OTTHUMP00000020967; terra-like protein
Protein Families:	Transcription Factors, Transmembrane



[View online »](#)

**Product images:**

TA334183 (2ug/ml) staining of paraffin embedded Human Kidney. Steamed antigen retrieval with citrate buffer pH 6, HRP-staining.