

Product datasheet for **TA334252**

FOXQ1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB, IHC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-FOXQ1 antibody: synthetic peptide directed towards the N terminal of human FOXQ1. Synthetic peptide located within the following region: MKLEVFVPRAAHGDKQGSDLEGAGGSDAPSPLSAAGDDSLGSDGDCAANS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	41 kDa
Gene Name:	forkhead box Q1
Database Link:	NP_150285 Entrez Gene 94234 Human Q9C009
Background:	FOXQ1 contains 1 fork-head DNA-binding domain. FOXQ1 mediates the interaction of Akt/protein kinase B with TGFb2. Foxq1 regulates differentiation of hair in satin mice.
Synonyms:	HFH1
Note:	Immunogen Sequence Homology: Human: 100%; Pig: 86%; Rat: 85%; Mouse: 85%
Protein Families:	Transcription Factors




[View online »](#)

Product images:



WB Suggested Anti-FOXQ1 Antibody Titration:
0.2-1 ug/ml; Positive Control: HepG2 cell lysate



Human Lung