

# **Product datasheet for TA334682**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### Parathyroid hormone related protein (PTHLH) Rabbit Polyclonal Antibody

#### **Product data:**

**Product Type:** Primary Antibodies

Applications:IHC, WBRecommended Dilution:WB, IHCReactivity:HumanHost:RabbitIsotype:IgG

Clonality: Polyclonal

**Immunogen:** The immunogen for anti-PTHLH antibody: synthetic peptide directed towards the middle

region of human PTHLH. Synthetic peptide located within the following region:

YKEQPLKTPGKKKKGKPGKRKEQEKKKRRTRSAWLDSGVTGSGLEGDHLS

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Purification: Protein A purified

Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

Predicted Protein Size: 20 kDa

**Gene Name:** parathyroid hormone-like hormone

Database Link: NP 002811

Entrez Gene 5744 Human

P12272





### Parathyroid hormone related protein (PTHLH) Rabbit Polyclonal Antibody - TA334682

**Background:** PTHLH is a member of the parathyroid hormone family. This hormone regulates

endochondral bone development and epithelial-mesenchymal interactions during the formation of the mammary glands and teeth. This hormone is involved in lactation possibly by regulating the mobilization and transfer of calcium to the milk. The receptor of this hormone, PTHR1, is responsible for most cases of humoral hypercalcemia of malignancy. Four alternatively spliced transcript variants encoding two distinct isoforms have been

observed.


**Synonyms:** BDE2; HHM; PLP; PTHR; PTHRP

**Note:** Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human:


100%; Rabbit: 100%; Guinea pig: 100%; Goat: 93%; Mouse: 93%; Bovine: 93%

**Protein Families:** Druggable Genome, Secreted Protein

## **Product images:**




Primary antibody dilution: 1:1000; Secondary antibody: Goat anti-rabbit HRP; Secondary antibody dilution:1:2000




WB Suggested Anti-PTHLH Antibody Titration: 2.0 ug/ml; Positive Control: HepG2 cell lysate





IHC Information:Rabbit Anti-PTHLH Antibody; Paraffin Embedded Tissue: Human Lung; Cellular Data: Alveolar cells; Antibody Concentration: 4.0-8.0 ug/ml; Magnification: 400X



Rabbit Anti-PTHLH Antibody; Paraffin Embedded Tissue: Human Kidney; Cellular Data: Epithelial cells of collecting tubule and renal tubule; Antibody Concentration: 4.0-8.0 ug/ml; Magnification: 400X