

Product datasheet for **TA334755**

ZNF179 (RNF112) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF179 antibody: synthetic peptide directed towards the middle region of human ZNF179. Synthetic peptide located within the following region: REFEEYVRQQDVATKRIFSALRVLPDTMRNLLSTQKDAILARHGVALLCK
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	68 kDa
Gene Name:	ring finger protein 112
Database Link:	NP_009079 Entrez Gene 7732 Human Q9ULX5
Background:	The specific function of this protein remains unknown. This gene encodes a member of the RING finger protein family of transcription factors. The protein is primarily expressed in brain. The gene is located within the Smith-Magenis syndrome region on chromosome 17.
Synonyms:	BFP; ZNF179
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Guinea pig: 100%; Mouse: 93%; Bovine: 93%; Rabbit: 93%



[View online »](#)

Protein Families: Transmembrane

Product images:



WB Suggested Anti-ZNF179 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:12500; Positive Control:
OVCAR-3 cell lysate