

Product datasheet for **TA334824**


ANKRD18A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ANKRD18A antibody is: synthetic peptide directed towards the C-terminal region of Human ANKRD18A. Synthetic peptide located within the following region: MRNVQEKCEKLEKDKKMLEEEVLNLKTHMEKDMVELGKLQEYKSELDERA
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	23 kDa
Gene Name:	ankyrin repeat domain 18A
Database Link:	NP_671728 Entrez Gene 253650 Human Q8IVF6
Background:	The function of this protein remains unknown.
Synonyms:	FLJ35740; FLJ40632; KIAA2015; MGC176674
Note:	Immunogen Sequence Homology: Human: 100%; Dog: 90%; Horse: 75%



[View online »](#)

Product images:


Host: Rabbit; Target Name: ANKRD18A; Sample
Tissue: Fetal Liver lysates; Antibody Dilution: 1.0
ug/ml