

## Product datasheet for **TA334830**


### CDRT15L2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-CDRT15L2 antibody is: synthetic peptide directed towards the N-terminal region of Human CDRT15L2. Synthetic peptide located within the following region: RLIPHPRRLWPFVRRRTQVPQDSPGQALAGQATPEIPSGPLPHIVLVQEE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	30 kDa
Gene Name:	CMT1A duplicated region transcript 15-like 2
Database Link:	<a href="#">NP_001177719</a> <a href="#">Entrez Gene 256223 Human</a> <a href="#">A8MXV6</a>
Background:	The function of this protein remains unknown.
Synonyms:	CDRT15L2
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transmembrane



[View online »](#)

**Product images:**

Host: Rabbit; Target Name: CDRT15L2; Sample  
Tissue: Jurkat Whole Cell lysates; Antibody  
Dilution: 1.0 ug/ml