

## Product datasheet for **TA334863**


### STAC2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-STAC2 antibody is: synthetic peptide directed towards the N-terminal region of Human STAC2. Synthetic peptide located within the following region: TEMSEKENEPDDAATHSPPGTVSALQETKLQRFKRSLSLKILRSKSLEN
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	44 kDa
Gene Name:	SH3 and cysteine rich domain 2
Database Link:	<a href="#">NP_945344</a> <a href="#">Entrez Gene 342667 Human</a> <a href="#">Q6ZMT1</a>
Background:	The function of this protein remains unknown.
Synonyms:	24b2; STAC2
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%



[View online »](#)

**Product images:**

Host: Rabbit; Target Name: STAC2; Sample  
Tissue: 293T Whole Cell lysates; Antibody  
Dilution: 1.0 ug/ml