

Product datasheet for **TA334892**


C3orf54 (FAM212A) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-FAM212A antibody is: synthetic peptide directed towards the middle region of Human FAM212A. Synthetic peptide located within the following region: RAPVASVPPVHHPRPKSTPDACLEHWQGLEAEDWTAALLNRGRSRQPLVL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	31 kDa
Gene Name:	family with sequence similarity 212 member A
Database Link:	NP_976248 Entrez Gene 389119 Human Q96EL1
Background:	The function of this protein remains unknown.
Synonyms:	C3orf54; INKA1
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Zebrafish: 85%



[View online »](#)

Product images:

Host: Rabbit; Target Name: FAM212A; Sample
Tissue: Fetal Liver lysates; Antibody Dilution: 1.0
ug/ml