

Product datasheet for **TA334911**

C3orf72 (FOXL2NB) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-C3orf72 antibody is: synthetic peptide directed towards the N-terminal region of Human C3orf72. Synthetic peptide located within the following region: PKMCLHMAVRHKAQKTGPGILQQRQKPPAPRASGGPALLGKRRGCSEAG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	18 kDa
Gene Name:	FOXL2 neighbor
Database Link:	NP_001035150 Entrez Gene 401089 Human Q6ZUU3
Background:	The function of this protein remains unknown.
Synonyms:	C3orf72
Note:	Immunogen Sequence Homology: Human: 100%



[View online »](#)

Product images:



Host: Rabbit; Target Name: C3orf72; Sample
Tissue: PANC1 Whole Cell lysates; Antibody
Dilution: 1.0 ug/ml