

Product datasheet for TA334939


DYNLT4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-TCTEX1D4 antibody: synthetic peptide directed towards the middle region of human TCTEX1D4. Synthetic peptide located within the following region: VCSVLGPRAGQGVHVSRALWDVARDGLASVSYTNTSLFAVATVHGLYC
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	23 kDa
Gene Name:	Tctex1 domain containing 4
Database Link:	NP_001013654 Entrez Gene 343521 Human Q5JR98
Background:	TCTEX1D4 (also known as TCTEX1D4) belongs to the dynein light chain Tctex-type family. The exact function of TCTEX1D4 remains unknown.
Synonyms:	TCTEX1D4
Note:	Immunogen Sequence Homology: Human: 100%; Rat: 79%


[View online »](#)

Product images:



WB Suggested Anti-TCTEX1D4 Antibody Titration:
0.2-1 ug/ml; Positive Control: 721_B cell lysate