


Product datasheet for **TA335116**

TMED8 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-TMED8 antibody: synthetic peptide directed towards the middle region of human TMED8. Synthetic peptide located within the following region: EIEEPVPAGDVERGSRSSLRGRYGEVMPVYRRDSHRDVQAGSHDYPGEGI
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	36 kDa
Gene Name:	transmembrane p24 trafficking protein family member 8
Database Link:	NP_998766 Entrez Gene 283578 Human Q6PL24
Background:	The specific function of the protein remains unknown.
Synonyms:	FAM15B
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%


[View online »](#)

Product images:


WB Suggested Anti-TMED8 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:1562500; Positive
Control: Jurkat cell lysate