

Product datasheet for **TA335231**

GPD2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-GPD2 antibody: synthetic peptide directed towards the N terminal of human GPD2. Synthetic peptide located within the following region: DILVIGGGATGSGCALDAVTRGLKTALVERDDFSSGTSSRSTKLIHGGVR
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	81 kDa
Gene Name:	glycerol-3-phosphate dehydrogenase 2
Database Link:	NP_000399 Entrez Gene 2820 Human P43304
Background:	The protein encoded by this gene localizes to the inner mitochondrial membrane and catalyzes the conversion of glycerol-3-phosphate to dihydroxyacetone phosphate, using FAD as a cofactor. Along with GDP1, the encoded protein constitutes the glycerol phosphate shuttle, which reoxidizes NADH formed during glycolysis. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jan 2010]
Synonyms:	GDH2; GPDM; mGPDH



[View online »](#)

Note: Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Yeast: 93%; Zebrafish: 93%

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Glycerophospholipid metabolism

Product images:



WB Suggested Anti-GPD2 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive Control: OVCAR-3 cell lysate