

Product datasheet for **TA335271**

PRRG1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-PRRG1 antibody: synthetic peptide directed towards the N terminal of human PRRG1. Synthetic peptide located within the following region: ECKEEFCTFEEAREAFENNEKTKEFWSTYTKAQQGESNRGSDWFQFYLT
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	25 kDa
Gene Name:	proline rich and Gla domain 1
Database Link:	NP_000941 Entrez Gene 5638 Human O14668
Background:	The exact function of PRRG1 remains unknown.
Synonyms:	PRGP1
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Rabbit: 100%; Rat: 93%; Mouse: 93%; Guinea pig: 93%
Protein Families:	Transmembrane



[View online »](#)

Product images:



WB Suggested Anti-PRRG1 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive Control: Human heart