

Product datasheet for **TA335323**


ANKRD18B Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ANKRD18B Antibody is: synthetic peptide directed towards the C-terminal region of Human ANKRD18B. Synthetic peptide located within the following region: QVNIPDIIHIHKEALTKVMESRQHVAEGKTEVQRLMMSESQEQDFFGHCG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	149 kDa
Gene Name:	ankyrin repeat domain 18B
Database Link:	NP_001231681 Entrez Gene 441459 Human
Background:	The function of this protein remains unknown.
Synonyms:	bA255A11.3; bA255A11.5



[View online »](#)

Product images:


Host: Rabbit; Target Name: ANKRD18B; Sample
Tissue: HepG2 Whole Cell lysates; Antibody
Dilution: 1.0 ug/ml