

Product datasheet for **TA335382**

DNAJC18 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-DNAJC18 Antibody is: synthetic peptide directed towards the N-terminal region of Human DNAJC18. Synthetic peptide located within the following region: EWTQTRQGEGNSTYSEEQLLGVQRIKKCRNYYEILGVSRDASDEELKKAY
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	41 kDa
Gene Name:	DnaJ heat shock protein family (Hsp40) member C18
Database Link:	NP_689899 Entrez Gene 202052 Human Q9H819
Background:	The function of this protein remains unknown.
Synonyms:	MGC29463
Note:	Immunogen Sequence Homology: Human: 100%; Bovine: 100%; Dog: 93%; Pig: 93%; Rat: 93%; Horse: 93%; Goat: 86%; Mouse: 86%; Rabbit: 86%; Guinea pig: 86%



[View online »](#)

Product images:



Host: Rabbit; Target Name: DNAJC18; Sample
Tissue: 293T Whole Cell lysates; Antibody
Dilution: 1.0 ug/ml