

Product datasheet for **TA335385**

C19orf25 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-C19orf25 Antibody: synthetic peptide directed towards the n terminal of human C19orf25. Synthetic peptide located within the following region: MGSKAKKRVLLPTRPAPPTVEQILEDVRGAPAEDPVFTILAPEDPPVPFR
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	18 kDa
Gene Name:	chromosome 19 open reading frame 25
Database Link:	NP_689695 Entrez Gene 148223 Human Q9UFG5
Background:	The specific function of the protein remains unknown.
Synonyms:	FLJ36666
Note:	Immunogen Sequence Homology: Pig: 100%; Human: 100%; Bovine: 100%; Rat: 93%; Mouse: 93%; Rabbit: 86%; Guinea pig: 86%



[View online »](#)

Product images:



WB Suggested Anti-C19orf25 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive
Control: OVCAR-3 cell lysate