

Product datasheet for **TA335398**

DCST2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-DCST2 Antibody is: synthetic peptide directed towards the C-terminal region of Human DCST2. Synthetic peptide located within the following region: AQRKDPEQAWLLQQQLQEVLGRSLSMESTSESSDLDEEKGPQQRKHGQQP
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	53 kDa
Gene Name:	DC-STAMP domain containing 2
Database Link:	NP_653223 Entrez Gene 127579 Human Q5T1A1
Background:	The function of this protein remains unknown.
Synonyms:	RP11-307C12.4
Note:	Immunogen Sequence Homology: Human: 100%; Rabbit: 77%
Protein Families:	Transmembrane



[View online »](#)

Product images:



Host: Rabbit
Target Name: DCST2
Sample Tissue: PANC1 Cell Lysate
Antibody Dilution: 1.0 µg/ml

Host: Rabbit; Target Name: DCST2; Sample
Tissue: PANC1 Whole cell lysates; Antibody
Dilution: 1.0 µg/ml