


Product datasheet for TA335439

REXO5 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-LOC81691 Antibody: synthetic peptide directed towards the middle region of human LOC81691. Synthetic peptide located within the following region: AEGGCCVMDELVKPENKILDYLTSGITKKILNPVTTKLKDVQRQLKAL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	85 kDa
Gene Name:	RNA exonuclease 5
Database Link:	NP_112203 Entrez Gene 81691 Human Q96IC2
Background:	The exact functions of LOC81691 remain unknown.
Synonyms:	DKFZp434J0315
Note:	Immunogen Sequence Homology: Human: 100%; Dog: 93%; Rat: 86%; Horse: 86%; Bovine: 86%; Rabbit: 86%; Mouse: 85%; Pig: 75%; Guinea pig: 75%


[View online »](#)

Product images:

WB Suggested Anti-LOC81691 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive Control: MCF7 cell lysate. LOC81691 is strongly supported by BioGPS gene expression data to be expressed in Human MCF7 cells