

Product datasheet for **TA335513**

TCPII Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-TCPII Antibody: synthetic peptide directed towards the middle region of human TCPII. Synthetic peptide located within the following region: DMVNYTIQSLQPHLQEHSIQYERAKFQELLNKQPSLLNHTTKWLTQAAGD
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	49 kDa
Gene Name:	t-complex 11
Database Link:	NP_061149 Entrez Gene 6954 Human Q8WWU5
Background:	TCPII may play an important role in sperm function and fertility.
Synonyms:	D6S230E; FPPR



Note: Immunogen Sequence Homology: Dog: 100%; Horse: 100%; Human: 100%; Pig: 93%; Rat: 93%; Mouse: 93%; Bovine: 93%; Rabbit: 93%; Guinea pig: 92%

Product images:



WB Suggested Anti-TCP11 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
Control: DU145 cell lysate