

Product datasheet for **TA335908**

Paox Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-Paox Antibody is: synthetic peptide directed towards the middle region of Rat Paox. Synthetic peptide located within the following region: FQLAAEFGLLGEKELSEENQLVETGGHVALPSVCTSSGTSVSLELVTEM
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	58 kDa
Gene Name:	polyamine oxidase (exo-N4-amino)
Database Link:	NP_001099781 Entrez Gene 293589 Rat
Background:	The function of this protein remains unknown.
Synonyms:	DKFZp434J245; MGC45464; PAO
Note:	Immunogen Sequence Homology: Human: 100%; Pig: 92%; Horse: 92%; Rat: 86%; Mouse: 86%; Bovine: 85%; Rabbit: 85%; Guinea pig: 85%



[View online »](#)

Product images:



Host: Rabbit; Target Name: Paox; Sample Tissue: Rat Pancreas lysates; Antibody Dilution: 1.0 ug/ml