

Product datasheet for **TA335991**


ZDHHC19 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ZDHHC19 Antibody: synthetic peptide directed towards the middle region of human ZDHHC19. Synthetic peptide located within the following region: AAGLLVPLSLLLLLIQALSVSADRTYK GKCRHLQGYNPFDQGCASNWYLT
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	26 kDa
Gene Name:	zinc finger DHHC-type containing 19
Database Link:	NP_001034706 Entrez Gene 131540 Human Q8WVZ1
Background:	ZDHHC19 belongs to the DHHC palmitoyltransferase family. It contains 1 DHHC-type zinc finger. The exact function of ZDHHC19 remains unknown.
Synonyms:	DHHC19
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transmembrane



[View online »](#)

Product images:

WB Suggested Anti-ZDHHC19 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
Control: Hela cell lysate