

Product datasheet for **TA336136**

LST 3TM12 (SLCO1B7) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-LST-3TM12 Antibody: synthetic peptide directed towards the middle region of human LST-3TM12. Synthetic peptide located within the following region: RAFFGLKVALIFPVLVLLTVFIFVVRKSKSHGKDTKVLENERQVMDEANLE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	71 kDa
Gene Name:	solute carrier organic anion transporter family member 1B7 (putative)
Database Link:	NP_001009562 Entrez Gene 338821 Human G3V0H7
Background:	The exact function of LST-3TM12 remains unknown.
Synonyms:	LST-3; LST-3TM12; LST3; SLC21A21
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transmembrane



[View online »](#)

Product images:



WB Suggested Anti-LST-3TM12 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 312500; Positive
Control: 721_B cell lysate