

Product datasheet for **TA336147**

C14orf180 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-C14orf180 Antibody: synthetic peptide directed towards the N terminal of human C14orf180. Synthetic peptide located within the following region: EDNRKCPPSILKRSRPEHHRPEAKPQRTSRRVWFREPPAVTVHYADKNA
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	18 kDa
Gene Name:	chromosome 14 open reading frame 180
Database Link:	NP_001008404 Entrez Gene 400258 Human Q8N912
Background:	C14orf180 is a multi-pass membrane protein. The function of the C14orf180 protein is not known.
Synonyms:	C14orf77; NRAC
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transmembrane



[View online »](#)

Product images:


WB Suggested Anti-C14orf180 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive
Control: Transfected 293T