

Product datasheet for **TA336204**

NKIRAS2 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-NKIRAS2 Antibody: synthetic peptide directed towards the middle region of human NKIRAS2. Synthetic peptide located within the following region: KKEVTIVVLGNKCDLQEQRVDPDVAQHWAKSEKVKLWEVSVADRRSLLE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	21 kDa
Gene Name:	NFKB inhibitor interacting Ras like 2
Database Link:	NP_001001349 Entrez Gene 28511 Human Q9NYR9
Background:	NPAL2 is a multi-pass membrane protein and it belongs to the NIPA family. The exact function of NPAL2 remains unknown. Western blots using two different antibodies against two unique regions of this protein target confirm the same apparent molecular weight in our tests.
Synonyms:	kappaB-Ras2; KBRAS2



[View online »](#)

Note: Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 86%

Product images:



WB Suggested Anti-NKIRAS2 Antibody Titration: 0.2-1 ug/ml; Positive Control: 721_B cell lysate. NKIRAS2 is strongly supported by BioGPS gene expression data to be expressed in Human 721_B cells